


ML2001

general purpose 3-in-1

pH



Order number ML2001
former order number -

Specifications
pH range 0 - 14
zero point (pH) 7.00 ± 0.25

temperature range 0 - 100 °C
temperature probe NTC 30 kΩ

shaft material glass
membrane glass HA
membrane resistance < 600 MΩ

reference system ARGENTHAL™, Ag⁺ trap
reference diaphragm ceramic
reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head MultiPin™
cable length -
connector -

cell constant (cm⁻¹) -
cell type -

indicating metal -

Storage 3 mol/L KCl

Special features
integrated NTC 30 kΩ temperature probe
ARGENTHAL™ reference system with Ag⁺ trap

Choice of cables
All Mettler-Toledo meters have BNC connectors.
Recommended cable: BNC/RCA (Cinch) 1.2 m (52300009)

Our cable selector guide offers a large choice of MultiPin™ cables for other connectors.

Scale 1:1