


ML2003

specialist 2-in-1

pH



former order number -

Specifications

pH range 1 - 11

zero point (pH) 7.00 \pm 0.25

temperature range 0 - 80 °C

temperature probe -

shaft material polysulfone

membrane glass LoT

membrane resistance < 150 M Ω

reference system Ag / AgCl

reference diaphragm porous PTFE

reference electrolyte gel

bridge electrolyte -

electrode head S7 (grey)

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage 3 mol/L KCl**Special features**

large porous PTFE reference junction

Choice of cables

All Mettler-Toledo meters have BNC connectors.

Recommended cable: BNC 1.2m (52300004)

Our cable selector guide offers a large choice of S7 cables for other connectors. Grey S7 connectors are compatible with red S7 connectors. The color is only design-related.

Scale 1:1